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**Menurut dokumen CLSI EP17-A2, Vol. 32, No.8

Rentang Rujukan [15]

Pedoman National Cholesterol Education Program (NCEP):

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Pustaka

- Grundy SM et al. *Guide on the Management of Blood Cholesterol: A Report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines*. Circulation. 2018; 138:e1082-e1143.
- Gordon DJ, Probstfield JL, Garrison RJ, et al. High density lipoprotein cholesterol and cardiovascular disease. Four prospective American Studies. Circulation. 1989;79: 8-15.
- Favari E, Chroni A, Tietge UJF et al. High Density Lipoproteins: From Biological Understanding to Clinical Exploitation. Springer Verlag; Volume 224, 2015; p. 161-206.
- Barter PJ, Nicholls S, Rye KA, Anantha-ramaiah G, Navab M and Fogelman AM (2004). Antiinflammatory properties of HDL. Circ. Res. 95.
- Lee JS, Chang P-Y, Zhang Y, Kizer JR, Best LG and Howard BV. Triglyceride and HDL-C Dyslipidemia and Risks of Coronary Heart Disease and Ischemic Stroke by Glycemic Dysregulation Status: The Strong Heart Study. Diabetes Care 2017; 40: 529-537.
- Chapman, M John et al. "Triglyceride-rich lipoproteins and high-density lipoprotein cholesterol in patients at high risk of cardiovascular disease: evidence and guidance for management". European heart journal volume 32, 11 (2011): 1345-61. 764-772.
- Rifai N, Warnick GR. Lipids, Lipoproteins, Apolipoproteins, and Other Cardiovascular Risk Factors, In: Burtis CA, Ashwood ER and Burns DE, editors. Tietz Textbook of Clinical Chemistry. 4th ed, Missouri. Elsevier Saunders company; 2006. p. 903-981.
- Recommendation of the Second Joint Task Force of European and other Societies on Coronary Prevention. Prevention of coronary heart disease in clinical practice. Eur Heart J 1998; 19: 1434-503.
- Langlois MR, Blaton VH. Historical milestones in measurement of HDLcholesterol: Impact on clinical and laboratory practice. Clin Chimica Acta 2006; 369: 168-178.
- Milda T, Nishimura K, Okamura T, et al. Validation of homogeneous assays for HDL-cholesterol using fresh samples from healthy and disease subjects. Atherosclerosis 2014;233(1); 253-9.
- Bakker AJ, Mucke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007; 45 (9): 1240-1243.
- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 22-3.
- Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Volume 1 and 2. Washington DC: The American Association for Clinical Chemistry Press 2000.
- Young DS. Effects on Clinical Laboratory Test – Drugs Disease, Herbs & Natural Products, accessed on May 2020. Published by AACC Press and John Wiley and Sons, Inc..
- Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults. Executive Summary of the Third Report of The National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (adult Treatment Panel III). JAMA 2001; 285(19) 2486-2497.

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