

# np nerbe plus

GERMANY

np GREEN TIP  
HANDLING. LIQUID.



## KOMPATIBILITÄTSTABELLE PIPETTENSPIZTEN

### COMPATIBILITY CHART OF PIPETTE TIPS

	10 µl short	10 µl extra long	0,1–20 µl super slim	100–200 µl
Pipettor Brand / Model	06-360-2008 06-360-2018	06-361-2008 06-361-2018	06-362-2008 06-362-2018 06-372-2018	06-365-2008 06-365-2018 06-375-2018
Sartorius Biohit M-Line™		M10	M10	M100, M200
Sartorius Biohit Proline™	0,1–3 µl, 0,5–10 µl	0,1–3 µl, 0,5–10 µl	0,1–3 µl, 0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Sartorius Biohit Proline Plus, Picus, Picus Nxt, Tacta™	0,1–3 µl, 0,5–10 µl, 0,2–10 µl	0,1–3 µl, 0,5–10 µl, 0,2–10 µl	0,1–3 µl, 0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl, 5–120 µl
Sartorius Biohit Proline Electronic	0,1–3 µl, 0,5–10 µl	0,1–3 µl, 0,5–10 µl	0,1–3 µl, 0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Brand Transferpette S™	0,5–10 µl	0,5–10 µl	0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Brand Transferpette Electronic™	0,5–10 µl	0,5–10 µl	0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Capp™ ecopipette / Aero(singlechannel)	0,2–2 µl, 0,5–10 µl	0,2–2 µl, 0,5–10 µl	0,2–2 µl, 0,5–10 µl	2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl
Capp™ electronic / AERO multichannel	0,5–10 µl, 2–20 µl	0,5–10 µl, 2–20 µl	0,5–10 µl, 2–20 µl	5–50 µl, 10–100 µl, 20–200 µl
CAPP Bravo	0,2–2 µl, 0,5–10 µl	0,2–2 µl, 0,5–10 µl	0,2–2 µl, 0,5–10 µl	2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl
Eppendorf Reference™		1–10 µl, 2–20 µl	1–10 µl, 2–20 µl	2–20 µl, 10–100 µl, 50–200 µl
Eppendorf Research™		0,1–2,5 µl, 1–10 µl	0,1–2,5 µl, 1–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Eppendorf Research Plus™	0,1–2,5 µl, 0,5–10 µl	0,1–2,5 µl, 0,5–10 µl	0,1–2,5 µl, 0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Eppendorf Explorer™ und Explorer Plus™		0,5–10 µl	0,5–10 µl	10–100µl
Finnpipette Electronic™	1–10 µl	1–10 µl	1–10 µl	5–50 µl, 20–200 µl
Finnpipette, Digital, F1 & F2™	1–10 µl	1–10 µl	1–10 µl	5–50 µl, 20–200 µl
Finnpipette Focus	0,3–3 µl, 0,5–5 µl, 1–10 µl	0,5–5 µl, 1–10 µl	0,5–5 µl, 1–10 µl	
Gilson Pipetman Classic / G / L / M™	P2, P2G, P2L, P10, P10G, P10L, P10M	P2, P2G, P2L, P10, P10G, P10L, P10M	P2, P2G, P2L, P10, P10G, P10L, P10M	P2, P2G, P2L, P20, P20G, P20L, P20M, P200, P200G, P200L, P200M
Gilson Pipetman Ultra™				U20, U200
Hamilton™	0,2–2 µl, 1–10 µl	0,2–2 µl, 1–10 µl	0,2–2 µl, 1–10 µl	2,5–25 µl, 10–100 µl, 30–300 µl
Labnet BioPette, BioPiette A, Labpette, E™	0,1–2 µl, 0,5–10 µl	0,1–2 µl, 0,5–10 µl	0,1–2 µl, 0,5–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Nichiryo Nichipet EX™	0,2–2 µl, 1–10 µl	0,2–2 µl, 1–10 µl	0,2–2 µl, 1–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Socorex Acura & Calibra 822™				2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl
Biohit Proline® Plus Single Channel Pipettes	0,1–3 µl, 1–10 µl	0,1–3 µl, 1–10 µl	0,1–3 µl, 1–10 µl	2–20 µl, 10–100 µl, 20–200 µl
Biohit eLine Advanced Electronic Pipette				5–120 µl
Labnet / HTL – Labpette Discovery Digital Pipette	0,2–2 ul	0,2–2 ul	0,2–2 ul	2–20 µl, 10–100 µl, 20–200 µl

	300 µl	1,000 µl	1,250 µl
Pipettor Brand / Model	06-367-2008 06-367-2018 06-377-2018	06-368-2008 06-368-2018	06-369-2008 06-369-2018 06-379-2018
Sartorius Biohit M-Line™	M20, M100, M200	M1,000	M1,000
Sartorius Biohit Proline™	2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl	100–1000 µl	100–1000 µl
Sartorius Biohit Proline Plus, Picus, Picus Nxt, Tacta™	2–20 µl, 10–100 µl, 20–200 µl	100–1000 µl	100–1000 µl
Sartorius Biohit Proline Electronic	5–1000 µl	10–500 µl, 50–1000 µl	
Brand Transferpette S™	20–200 µl, 15–300 µl	100–1,000 µl	100–1,000 µl
Brand Transferpette Electronic™	15–300 µl		
Capp™ ecopipette / Aero(singlechannel)	2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl	100–1,000 µl	100–1,000 µl
Capp™ electronic / AERO multichannel	5–50 µl, 10–100 µl, 20–200µl	100–1,200 µl	100–1,200 µl
CAPP Bravo	2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl	100–1,000 µl	100–1,000 µl
Eppendorf Reference™	10–100 µl, 50–200 µl	100–1,000 µl	100–1,000 µl
Eppendorf Research™	2–20 µl, 10–100 µl, 20–200 µl, 30–300 µl	100–1,000 µl	100–1,000 µl
Eppendorf Research Plus™	2–20 µl, 10–100 µl, 20–200 µl, 30–300 µl	100–1,000 µl	100–1,000 µl
Eppendorf Explorer™ und Explorer Plus™	2–20 µl, 10–100 µl, 20–200 µl, 30–300 µl	100–1,000 µl	100–1,000 µl
Finnpipette Electronic™	5–50 µl, 20–200 µl, 30–300 µl	100–1,000 µl	100–1,000 µl
Finnpipette, Digital, F1 & F2™	5–50 µl, 20–200 µl, 30–300 µl	100–1,000 µl, 200–1,000 µl	100–1,000 µl, 200–1,000 µl
Finnpipette Focus	5–50 µl, 30–300 µl	100–1,200 µl	100–1,200 µl
Gilson Pipetman Classic / G / L / M™	P20, P20G, P20L, P20M, P100, P100G, P100L, P100M, P200, P200G, P200L, P200M	P1,000, P1,000G, P1,000L, P1,000M	P1,000, P1,000G, P1,000L, P1,000M
Gilson Pipetman Ultra™	U20, U100, U200		
Hamilton™	10–100 µl, 30–300 µl	100–1,000 µl	100–1,000 µl
Labnet BioPette, BioPiette A, Labpette, E™	2–20 µl, 10–100 µl, 20–200 µl	100–1,200 µl	100–1,200 µl
Nichiryo Nichipet EX™	2–20 µl, 10–100 µl, 20–200 µl	100–1,000 µl	100–1,000 µl
Socorex Acura & Calibra 822™	2–20 µl, 5–50 µl, 10–100 µl, 20–200 µl	100–1,000 µl	100–1,000 µl
Biohit Proline® Plus Single Channel Pipettes	2–20 µl, 10–100 µl, 20–200 µl	100–1,000 µl	100–1,000 µl
Biohit eLine Advanced Electronic Pipette	5–120 µl		
Labnet / HTL – Labpette Discovery Digital Pipette	2–20 µl, 10–100 µl, 20–200 µl		

## OBERFLÄCHENOPTIMIERTE FILTERSPITZEN PREMIUM SURFACE FILTER TIPS

CE IVD

• Konformitätsgeprüft nach ISO 8655  
Conformity tested to ISO 8655




06-642-0000  
06-642-5300

06-662-0000  
06-662-5300

06-672-0000  
06-672-5300

REF	MAT							
06-642-0000	PP	53,4 mm	100 µl		1.000		10.000	
06-642-5300	PP	53,4 mm	100 µl	96		960	5.760	
06-662-0000	PP	53,4 mm	200 µl		1.000		10.000	
06-662-5300	PP	53,4 mm	200 µl	96		960	5.760	
06-672-0000	PP	59,35 mm	300 µl		1.000		10.000	
06-672-5300	PP	59,35 mm	300 µl	96		960	5.760	



06-693-0000  
06-693-5300

06-695-0000  
06-695-5300

REF	MAT							
06-693-0000	PP	83,5 mm	1.000 µl		500		5.000	
06-693-5300	PP	83,5 mm	1.000 µl	96		576	3.456	
06-695-0000	PP	98,2 mm	1.250 µl		500		5.000	
06-695-5300	PP	98,2 mm	1.250 µl	96		576	3.456	



nerbe plus GmbH & Co. KG  
Porschestraße 25  
D-21423 Winsen / Luhe  
Germany  
Tel: +49 (0) 4171 7849 0  
Fax: +49 (0) 4171 645 45  
E-Mail: [info@nerbe-plus.de](mailto:info@nerbe-plus.de)  
Web: [www.nerbe-plus.de](http://www.nerbe-plus.de)

