


MIKROSKOP BINOKULER BM-001

KEMENKES RI AKD 10203320225



Spesifikasi BM-001 1600x :

Lensa Okuler	: Binokuler, Kemiringan 30°, Jarak interpupillary 55-75 mm, WF 10x18 mm
Revolver	: 4 Posisi Lensa Obyektif (Quadruple)
Lensa Obyektif	: PLAN Objective 4X,10X,40X,100X (oil)
Meja Kerja	: 2 Tingkat Mekanikal ukuran 140x135 mm / 75x50 mm
Kondensor	: Abbe NA 1.25 dengan Iris Diafragma dan Filter
Sistem Fokus	: Sistem pengaturan makrometer dan mikrometer, 37.7 mm per rotasi, fine stroke 0.2 mm per rotasi, fine division 0.002 mm, moving range 28 mm, safety autofocus stop unit
Sumber Cahaya	: Lampu LED 3W, tingkat pencahayaan dapat disesuaikan, plan-concave mirror
Lensa Filter	: Biru dan Hijau
Kemasan	: Sterofoam, Dus
Dimensi	: P x L x T (40 x 30 x 50 cm)
Berat Netto	: 8 Kg
Berat Gross	: 12 Kg
Aksesoris	: WF16x/12 eyepiece (23mm Ø) Cleaning kit Mikroskop N-PLAN objective 20x/0.40 USB Camera Eyepiece C-E2 ADVAN TBOOK Intel N100 Ram 4GB 128GB 14" HD Windows 11



Didistribusikan oleh :

PT. GERMARINDO BIOTEKNIKA

Jl. Ciledug Raya No. 232 Lt.3, Jakarta Selatan - DKI. Jakarta, Indonesia

 germarindobioteknika@gmail.com  +6221 722 7169