

BINOCULAR MICROSCOPE

BM-001



Specification BM-001 :

Viewing Head	: Compensation free binocular head inclined at 30° Inclined Interpupillary (55-75mm)
Eyepiece	: Wide Field eyepiece WF10X18mm
Nosepiece	: Quadruple Nosepiece
Objective	: PLAN Objective 4X, 10X, 40X, 100X (oil)
Stage	: Double Layers Mechanical Stage 140X135mm/ 75X50mm
Condenser	: Abbe NA1.25 with Iris Diaphragm & Filter
Focusing System	: Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 28mm, Safety autofocus stop unit.
Illumination	: 3W LED Illumination, Brightness Adjustable Plan-concave Mirror
Package	: Styrofoam Carton Wooden Case (optional)
Filter	: Blue & Green Filter
Package Size	: LxWxH (40x30x50cm)
Nett Weight	: 8 kg
Gross Weight	: 12 kg